


B.Sc.-2ND Year

CHEMISTRY - 3RD Semester

INDEX

S.r. No.	Name of the Experiment	Page No.	Date of Experiment	Date of Submission	Remarks
1	To determine the boiling point of the given ^{organic} compound.				
2	To determine the boiling point of the given organic compound.				
3	To determine the boiling point of the given organic compound.				
4	To determine the boiling point of the given organic compound.				
5	To determine the boiling point of the given organic compound.				
6	To determine the viscosity supplied solution using ostwald viscometer.				
7	To determine the surface tension of supplied solution using stalagmometer.				
8	To prepare 250 ml dilute solution of 0.001 M sulphuric acid. Provided 0.1 M sulphuric acid				
9.	To prepare 250 ml dil. Solution of 0.001 M H_2SO_4 . Provided 0.1 M H_2SO_4 .				

EX-1



Determination of Boiling Point

Experiment no. -1

Object -

To determine the boiling point of the given organic compound.

Apparatus -

Boiling tube, spirit lamp, thermometer, burette stand etc.

Procedure -

The liquid whose boiling point is to be determined is taken in a test tube or boiling tube.

The tube is clamped in a stand and heated gently. A thermometer is inserted in the tube so that its bulb remains at a distance of about 1 cm above the surface of the liquid. When the liquid begins to boil the temperature is noted. This is boiling point.

Result -

The given organic compound is boiling point -


Cyclohexanol - 161.8°C

Fixed value - 161.8°C

Observed value - 158.9°C

Teacher's Signature :

Ex-2



Determination of Boiling Point

Experiment no. -2

Object -

To determine the boiling point of the given organic compound.

Apparatus -

Boiling tube, spirit lamp, thermometer, burette stand etc.

Procedure -

The liquid whose boiling point is to be determined is taken in a test tube or boiling tube. The tube is clamped in a stand and heated gently. A thermometer is inserted at a distance of about 1 cm above the surface of the liquid. When the liquid begins to boil, the temperature is noted. This is boiling point.

Result -

The given organic compound is boiling point.


Cyclohexanone - 155.6°C

fixed boiling point - 155.6°C

Observed boiling point - 154.0°C

Teacher's Signature :

Ex-3



Determination of boiling point

Experiment no. -3

Object -

To determine the boiling point of the given organic compound.

Apparatus -

Boiling tube, spirit lamp, thermometer, burette stand etc.

Procedure -

The liquid whose boiling point is to be determined is taken in a test-tube or boiling tube. The tube is clamped in a stand and heated gently. A thermometer is inserted in the tube so that its bulb remains at a distance of about 1cm above the surface of the liquid. When the liquid begins to boil, temperature is noted. This is boiling point.

Result -

The given organic compound is boiling point -


Benzaldehyde - 178.1°C

Fixed m.p. - 178.1°C

Observed m.p. - 177.0°C

Teacher's Signature :

Ex-4



Determination of boiling point

Experiment no. - 4

Object -

To determine the boiling point of the given organic compound.

Apparatus -

Boiling tube, spirit lamp, thermometer, burette stand etc.

Procedure -

The liquid whose boiling point is to be determined is taken in a test tube or boiling tube. The tube is clamped in a stand and heated gently. A thermometer is inserted in the tube so that its bulb remains at a distance of about 1 cm above the surface of the liquid. When the liquid begins to boil, the temperature is noted. This is boiling point.

Result -

The given organic compound is boiling point -


(Acetylacetone) - 140°C

Fixed B.P. = 140°C

Observed B.P. = 138.8°C

Teacher's Signature :

Ex-5



Determination of boiling point

Experiment no- 5

Object -

To determine the boiling point of the given organic compound.

Apparatus -

Boiling tube, spirit lamp, thermometer, burette stand etc.

Procedure -

The liquid whose boiling point is to be determined is taken in a test tube or boiling tube. The tube is clamped in a stand and heated gently. A thermometer is inserted in the tube so that its bulb remains at a distance of about 1 cm above the surface of the liquid. When the liquid begins to boil, the temperature is noted. This is boiling point.

Result -

The given organic compound is boiling point-


108°C (Isobutyl Alcohol)

Fixed boiling point - 108°C

Observed boiling point - 107°C

Teacher's Signature :

Ex - 6



Ostwald viscometer

Experiment no-6

Object -

To determine the viscosity supplied solution using ostwald viscometer.

Apparatus -

Viscometer, beaker, glass rod, stop watch, chemical balance etc.

Procedure -

The ostwald viscometer is wash through with chromic acid rinsed with a large volume of water and finally with acetone (CH_3COCH_3) repeat the processes three or four times to insure that there is no serious deviation in any case from the average time. find out the mean time.

Observation table -

(a) Room temperature = 20°C

(b) Measurement of flow time

S.No.	Water	Supplied solution
	flow time, t_w	flow time, t_s
1.		
2	12	15
3	13	16
4	11	13
Mean value.	$\frac{12+13+11}{3} = 12$	$\frac{15+16+13}{3} = 14.66$

Teacher's Signature :

Calculation -

By substituting the values of t_1 , t_2 , d_1 , and n_w in the viscosity of supplied solution can be calculated.

$$\eta_{\text{solution}} = \frac{d_{\text{solution}} \times t_{\text{solution}}}{d_{\text{water}} \times t_{\text{water}}} \times \eta_{\text{water}}$$

$$= \frac{1.22 \times 14.66}{1 \times 12} \times 0.0101$$

$$\boxed{\eta_{\text{solution}} = 0.01505 \text{ Poise}}$$

Expt. No.....

Page No.....

Date...../...../.....

- (c) Density of water = $d_2 = 1 \text{ gm}$
(d) Viscosity of water = $0.0101 \text{ poise at } 20^\circ\text{C}$
(e) Density of benzene = d_1
(i) Weight of the empty pyknometer = $a \text{ g} = 19.18 \text{ gm}$
(ii) Weight of pyknometer filled with water = $b \text{ g} = 23.58 \text{ gm}$
(iii) Weight of water filled = $(b-a) \text{ gm} = 4.40 \text{ gm}$
(iv) Weight of pyknometer filled with supplied solution =
 $c \text{ g} = 24.58 \text{ gm}$
(v) Weight of solution filled = 5.40 gm

$$\text{Density of benzene} = \frac{\text{Weight of solution}}{\text{Weight of water}}$$


$$d_1 = \frac{(c-a)}{(b-a)} = \frac{5.4}{4.4} = 1.22 \text{ gm/ml}$$

Result-

Viscosity of supplied solution = 0.01505 poise

Teacher's Signature :

Ex - 7



Stalagmometer

Experiment no. 7

Object -

To determine the surface tension of supplied solution using stalagmometer.

Apparatus -

Stalagmometer, pyknometer or specific gravity bottle, supplied solution, water, beaker, thermometer.

Procedure -

Clean the stalagmometer, thoroughly with chromic acid and then rinse it with a large volume of water several times. Dry the stalagmometer. Dip one end of the stalagmometer in a beaker containing water. Such the water till its level rises 1-2 cm above the mark A, keep the stalagmometer vertical by fixing it in a clamp and allow the drops to fall and start counting the number of drops till the water level reaches the lower mark A. Let the number of drops be n_w . Now dry the stalagmometer either by blowing hot air or in an electric oven, and repeat the above procedure with benzene. Let the number of drops be n_i .

Teacher's Signature :

Calculation

With the help of the following equation calculate the surface tension of given solution,

$$\gamma_i = \frac{n_w}{n_i} \times \frac{d_i}{dw} \times \gamma_w$$

Substitute the value of γ_w , n_w , d_i , n_i and dw in the above equation and calculate the surface tension of benzene.

$$\gamma_i = \frac{88.67 \times 1.22 \times 72.8}{101 \times 1}$$

$$\gamma_i = 77.97 \text{ Dyne / cm}$$

Observation table-

(A) Data -

Temperature = 20°C Density of water, $d_w = 1\text{ gm/ml}$ at 20°C Surface tension of water, $\gamma_w = 72.8\text{ dynes/cm}$ at 20°C

(B) Table for the drop Numbers -

S.No.	No. of drops water	Given solution No. of drops
1.	81	85
2.	86	85
3.	86	86
	Mean(n_w) = 84.67	Mean(n_i) = 85.33

(c) Determination of density of given solution -

Weight of the empty pyknometer = $a\text{ g} = 18.18\text{ gm}$ Weight of pyknometer filled with water = $b\text{ g} = 22.5\text{ gm}$ $\therefore$ Weight of water = $(b-a)\text{ g} = 4.40\text{ gm}$ Weight of pyknometer filled with solution = $c\text{ g} = 23.58\text{ gm}$ $\therefore$ Weight of solution = $(c-a)\text{ gm} = 5.40\text{ gm}$

Density of solution = $\frac{\text{Weight of solution}}{\text{Weight of water}}$

$$= \frac{(c-a)}{(b-a)} d_w = \frac{5.40}{4.40} \times 1.22\text{ gm/ml}$$

Teacher's Signature :

Experiment no. 8

Object -

To prepare 250 ml dilute solution of 0.001 M sulphuric acid. Provided 0.1 M sulphuric acid.

Theory -

For preparing dilute solution of the required concentration, the molarity equation $m_1 V_1 = m_2 V_2$ is used.

To prepare 250 ml of 0.001 M sulphuric acid the volume of 0.1 M sulphuric acid required will be calculated as follows -

$$0.1 \text{ M} \times V_1 = 0.001 \text{ M} \times 250$$

$$V_1 = \frac{0.001 \times 250}{0.1} = 2.5 \text{ ml}$$

Thus, 2.5 ml of 0.1 M sulphuric acid is taken and diluted to 250 ml in a measuring flask.

Apparatus -

Measuring flask, graduated pipette, measuring cylinder, beakers, glass rod, 0.1 M sulphuric acid.

Procedure -

[i] The volume of 0.1 M H_2SO_4 required for preparing 250 ml of 0.001 M concentration is 2.5 ml.

[ii] Take a 250 ml measuring flask and put 100 ml of distilled water in it.

Teacher's Signature :

- [iii] Transfer with the help of a graduate pipette 2.5 ml of 0.1 M H_2SO_4 into the measuring flask slowly with stirring.
- [iv] After adding all the acid add more of distilled water to make the total solution to 250 ml.

A Mole —

Chemists have defined a mole as a 'definite number' of particles (atoms, molecules, ions or electrons)

This definite number is called Avogadro number, equal to 6.022×10^{23} , in honour of Amadeo Avogadro.

At present a mole is defined as the amount of a substance containing as many atoms, molecule, ions, electrons or other elementary quantities as there are carbon atoms in exactly 12 gm of ^{12}C .

The following are the definitions of 'mole' represented in the form of equation.

- [i] The number of moles of atoms — $\frac{\text{Weight in gm}}{\text{Atomic weight}}$
- [ii] Number of moles of molecule — $\frac{\text{weight in gm}}{\text{molecular weight}}$
- [iii] Number of moles of gasses — $\frac{\text{Volume at NTP}}{\text{Standard molar volume}}$

(Standard molar volume is the volume occupied by 1 mole of any gas at NTP which is equal to 22.4 liters.)

Teacher's Signature :

[iv]
$$\text{Number of moles of atoms/molecules/ions/electrons} = \frac{\text{No. of atom/molecule/ions/electrons}}{\text{Avagadro no.}}$$

[v]
$$\text{Number of moles of solution} = \text{molarity} \times \text{volume of solution in liters}$$

or
$$\text{No. of millimoles} = \text{Molarity} \times \text{Volume in ml}$$

$$= \frac{\text{Millimoles}}{1000} = \text{moles}$$

[vi] for a compound M_x , x moles of $N = y$ moles of M

Molecular weight -

Molecular weight of a compound is defined as the weight of a molecules of a compound relative to a carbon atom, the atomic weight of which is supposed to be exactly 12. The molecular weight when expressed in gram is called gram molecular weight. The molecular weight in gram is the weight of 1 mole of molecules, example -

Example - Molecular weight of O_2 is 32 gm is the weight of 1 mole oxygen molecules.

$$\text{Molecular weight} = \frac{\text{Weight of molecules in gram}}{\text{No. of moles of molecules}}$$

Teacher's Signature :

Experiment no. - 9

Object -

To prepare 250 ml dilute solution of 0.001 M sulphuric acid. Provided 0.1 M sulphuric acid.

Theory -

For preparing dilute solution of the required concentration, the molarity equation $M_1V_1 = M_2V_2$ is used.

To prepare 250 ml of 0.001 M sulphuric acid the volume of 0.1 M sulphuric acid required will be calculate as follows —

$$0.1 \text{ M} \times V_1 = 0.001 \text{ M} \times 250$$

$$V_1 = \frac{0.001 \times 250}{0.1} = 2.5 \text{ ml}$$

Thus, 2.5 ml of 0.1 M sulphuric acid is taken and dilute to 250 ml in a measuring flask.

Apparatus —

Measuring flask, graduate pipette, measuring cylinder, beakers, glass rod, 0.1 M sulphuric acid.

Procedure — (i) The volume of 0.1 M H_2SO_4 required for preparing 250 ml of 0.001 M Concentration is 2.5 ml.

(ii) Take a 250 ml measuring flask and put 100 ml of distilled water in it.

Teacher's Signature :

(iii)

Transfer with the help of a graduated pipette 2.5 ml of 0.1 M H_2SO_4 into the measuring flask slowly with stirring.

(iv)

After adding all the acid add more of distilled water to make the total solution to 250 ml.

Normality -

Number of the gram-equivalent of a substance dissolved in 1000 cc and 1 liter of the solution is called the normality of the solution.

If one gram equivalent weight of the substance is dissolved in 1 litre of the solution, then its normality is equal to one and it is called normal solution. Then Normality is represented by N.

Thus,
$$N = \frac{\text{No. of gram equivalent of the solute}}{\text{Volume of the solution in litre.}}$$

or
$$N = \frac{\text{Weight in gram of the solute per litre of the solution}}{\text{Gram equivalent mass of the solute.}}$$

We can prepare the solution of desired normality by using the following relation.

$$W = \frac{E \cdot N \cdot V}{1000}$$

Where
 W = amount of the substance (Solute)
 E = Equivalent weight of the substance (Solute)
 V = Volume of the solution
 N = Normality of the solution.

Teacher's Signature :

Molality -

Number of moles of a substance dissolved in 1000 gm or 1 kg of the solvent is called the molality of the solution. If one mole of the substance is dissolved in 1 kg of the solvent, then its molality is equal to one and it is called molal solution.

Molality is represented by m .

Then,

$$m = \frac{\text{No. of moles of solute}}{\text{Mass of the solvent in kg}}$$

or

$$m = \frac{\text{Weight in grams of the solute per kg of the solvent}}{\text{Gram molecular mass of the solute.}}$$

Molarity -

No. of moles of a substance dissolved in one litre of the solution is called the molarity of the solution. If the one litre of the solution contains one mole of the solute, then it is called molar solution. Molarity of the solution is represented by M .

Thus,

$$M = \frac{\text{No. of moles of the solute}}{\text{Volume of the solution in litre}}$$

or

$$M = \frac{\text{Weight in grams of the solute per litre of the solution}}{\text{Gram molecular mass of the solute.}}$$

Experiment no-10

Object -

To calibrate the given burette.

Procedure -

1. Before calibrate, it is necessary that the burette is thoroughly cleaned. Thus, the burette is first filled with cleaning mixture which consist of concentrated sulphuric acid and potassium-dichromate and allowed to stand for half an hour. After that the cleaning mixture is poured off and thoroughly rinsed three to four times with tap water and finally with distilled water.
2. The burette is fixed vertically in a stand and distilled water is filled slightly above the zero mark. The air bubbles at the tip of burette is continuously removed by opening the stop clock to its full extent. Now the level of water in burette is brought to the zero mark by opening the stop clock. Wait for one minute to avoid drainage error. Then the stop clock is slowly turned so that liquid flows out dropwise into a weighted stoppered 100 ml Erlenmeyer flask. A period of 30 seconds should be required to deliver 5 ml. The level of water is brought as closely as possible to 5 ml mark. The tip of the burette is touched to the inside wall of the flask to remove adhering droplet of water and the stoppered flask is then weighed accurately. From the weight of water delivered, the volume delivered is

Teacher's Signature :

calculated and the result is compared with the observed volume as obtained by subtracting the first reading from the last. The difference between the first reading and the observed and calculating values gives the errors of the burette at the 5 ml graduation.

3. The burette is then filled to the zero mark again and water is allowed to flow out drop wise as before to the 10 ml marks at the rate not deliverer exceeding 10 ml per minute. The water is delivered in the same flask weight accurately. The difference of this mass to last mass gives the mass of next portion of water. Note the reading after one minute.

Observation -

Room temperature = 22°C

mass of empty Erlenmeyer flask = 10.65 gm

density of distilled water,


at 22°C = 0.9997 gm/ml

Obs. No.	Volume of water delivered (ml)	mass of water delivered (b)	mass of flask	mass of water in flask (a)	Error = (a-b)
1	50	2	10.65	12.65	2
2	50	3	10.00	13.00	3
3	50	2	11.00	13.00	2
4	50	2	10.65	12.65	2

Result - Error in the burette was = 2.25 ml

Teacher's Signature :

Ex-11



Apparatus for determining C.S.T

Experiment no-11

Object -

Determine the mutual solubility curve of phenol and water system and hence the critical solution temperature.

Theory -

In order to determine the mutual solubility curve of this system synthetic method is followed. known amounts of the two liquid are heated slowly, and a temperature is recorded at which the two layers just disappears. The temperature is once more determined at the appearance of turbidity when the system is cooled. The mean of the two temperature gives the miscibility temperature of the two liquids in that proportion.

Requirements -

Bailing tubes, stirrer, a large tube as an air jacket, thermometer reading to 0.1°C , pure phenol.

Observation -

Weight of phenol - 8.0 gm

Density of water - 1 gm/ml

Teacher's Signature :

Calculation -

A curve is plotted between miscibility temperature as ordinate and percentage by weight of phenol as abscissa. This curve is shown in the following figure and is known as mutual solubility curve of phenol water system. The maxima of the solubility curve will give the critical solution temperature.




fig. - Miscibility Curve